

## Accucold Freezer FS24LZ



1.4 cu. ft. capacity; temperature capability of -20°C

## Freezers, Standard Range: Accucold Freezer FS24LZ

---

Rating: Not Rated Yet

### Price

Sales price \$175.00

[Ask a question about this product](#)

### Description

The Accucold FS24LZ is a compact, general-purpose all-freezer designed for reliable cold storage in medical, scientific, and laboratory settings. With a 1.4 cu. ft. capacity and a temperature capability of -20°C, it is ideal for storing sensitive materials. This manual defrost freezer features a factory-installed lock, a reversible door, and a flat door liner compliant with CDC guidelines, ensuring secure and efficient operation. Its countertop-friendly size, seamless interior liner, and removable wire shelf enhance usability and ease of cleaning.

### Specifications:

Model: FS24LZ

Capacity: 1.4 cu. ft. (40 L)

Temperature Range: Capable of -20°C

Defrost Type: Manual

Convection Type: Static cooling

Door: Reversible with factory-installed lock

Door Liner: Flat, CDC-compliant

Shelving: 1 removable wire shelf

Interior Liner: Seamless for easy cleanup

Thermostat: Dial, located inside

Refrigerant: R600a (eco-friendly hydrocarbon-based)

Exterior Dimensions (H x W x D): 20.0" x 18.5" x 17.63" (51 cm x 47 cm x 45 cm)

Interior Dimensions (H x W x D): 14.38" x 13.5" x 13.25" (37 cm x 34 cm x 34 cm)

Depth with Door at 90°: 34.5" (88 cm)

Voltage/Frequency: 115 V AC/60 Hz

Amps: 0.8

Weight: 35 lbs (16 kg)

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)