

Across International Jacketed Glass Reactor R50





With Thermo Insulation Jacket, 50 L, 110V, 60 Hz, 140W.

Rating: Not Rated Yet

Price

Sales price \$3,750.00

[Ask a question about this product](#)

Description

Across International Jacketed Glass Reactor R50

Power :110V, 60 Hz, 1.7A, 140W

50L capacity

Specifications

Model: R50

Electrical requirements: 110V 60Hz, 140W

Reaction glass vessel capacity: 13 gallons (50L)

Cooling/heating jacket capacity: 4.2 gallons (16L)

Free high quality thermos insulation jacket with open/close viewing window

Drain port ground clearance: 12.5 inches. Unit height (fully assembled): 88 inches

Dimensions: 63" H, 20" L, 35" W

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.

226 Miller Ave.

South San Francisco, CA 94080

650-871-5707 office

650-276-7487 fax

sales@biotechequipmentsales.com