

Agilent 6890N



Agilent 6890N Gas Chromatograph w/TCD and FID detectors, split/splitless and packed column injectors

Rating: Not Rated Yet

Price

Sales price \$17,500.00

[Ask a question about this product](#)

Manufacturer [Agilent](#)

Description

!!!SOLD!!!

This Agilent 6890N Network Gas Chromatograph w/TCD and FID detectors, split/splitless and packed column injectors has barely been used. It comes w/a loaded laptop containing Chemstation software.

The 6890 represents the culmination of all the engineering and technological advancement that was achieved from every GC product line up through the 5890. The integration of a new technology for Electronic Pressure Control (EPC) was introduced, along with EPC options for auxiliary gases if needed. Revolutionary new inlet and detector designs deliver performance and productivity enhancements over previous capabilities, while allowing for simplified, user-installable "swapable" detector system

Electronic Pressure Control on detector and injector, ambient temperature +4 °C to +450 °C, temperature set point resolution: 1 °C, supports 6 oven temps with 7 programs, Negative temps are allowed, LAN, RS232 and remote communication, two analog output channels of 1 mV, 1 V, and 10 V, for a total of 6 outputs are standard.

Please call us at (650) 871-5707 to check availability of this product!