


## Agilent Bio SEC-5 Guard Column - New in Box



Particle size: 5  $\mu\text{m}$ ; Pore size: 300  $\text{\AA}$

Rating: Not Rated Yet

**Price**

Sales price \$250.00

[Ask a question about this product](#)

Manufacturer [Agilent](#)

### Description

The Agilent Bio SEC-5 Guard column is a size-exclusion chromatography (SEC) guard column designed to protect and extend the lifetime of a main analytical SEC column in biomolecule applications. It uses 5  $\mu\text{m}$  spherical, fully-porous silica particles coated with a neutral hydrophilic polymeric layer, enabling stable and reproducible separations of proteins, peptides, and other macromolecules.

Specifications

Particle size: 5 µm

Pore size: 300 Å

Column dimensions: 4.6 mm × 50 mm length

pH operating range: approx. pH 2.0 to 8.5.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)