

Agilent Bioanalyzer 2100 model #G2938A



With chip vortexer and priming station. With PC and loaded software.

Rating: Not Rated Yet

Price

Sales price \$8,500.00

[Ask a question about this product](#)

Manufacturer [Agilent](#)

Description

Model #: G2938A

The 2100 Bioanalyzer system is an established automated electrophoresis tool for the sample quality control of biomolecules. The 2100 Bioanalyzer instrument, together with the 2100 Expert Software and Bioanalyzer assays, provide highly precise analytical evaluation of various samples types in many workflows, including next generation sequencing (NGS), gene expression, biopharmaceutical, and gene editing

research.

Digital data is provided in a timely manner and delivers objective assessment of sizing, quantitation, integrity and purity from DNA, RNA, and proteins. Minimal sample volumes are required for an accurate result, and the data may be exported in a many different formats for ease-of-use.

Specifications:

Ambient Operating Temperature: 5-40C

Dimensions (WxDxH): 162 mm x 412 mm x 290 mm

Electrical req: 50-60 Hz / 100-240 VAC

Power Consumption: 60 VA

Weight: 10 kg

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com