

Alcatel Pascal 1005 L



Alcatel Pascal 1005 SD direct-drive vacuum pump

Rating: Not Rated Yet

Price

Sales price \$1,025.00

[Ask a question about this product](#)

Manufacturer [Alcatel](#)

Description

This Alcatel Pascal model 1005 SD vacuum pump is direct-drive & single stage. It has a pumping speed of 3.8 cfm, and an ultimate total pressure of 5×10^{-2} mbar w/a closed gas ballast. This pump was used for only a single type of application and has been thoroughly tested by our technicians.

Biotech Equipment Sales, Inc. based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

[Biotech Equipment Sales, Inc.](#)
265 Miller Ave
South San Francisco, CA 94080
650-871-5707 office
650-276-7437 fax
sales@biotechequipmentsales.com