

Amersham AKTA Prime



Low pressure-1MPa max with flow rate of 0.1-50mL minimum chromatographic system, designed primarily for protein purification, the system is available for sale with a computer with PrimeView software installed for 9,800 dollars or without a computer and so

Rating: Not Rated Yet

Price

Sales price \$9,800.00

[Ask a question about this product](#)

Manufacturer [Amersham](#)

Description

This AKTA Prime has low pressure-1MPa max with flow rate of 0.1-50mL minimum chromatographic system, designed primarily for protein purification, the system is available for sale with a computer with **PrimeView** software installed for 9,800 dollars or without a computer and software but with a REC 112 chart recorder for 8,500 dollars.