

## Amersham Ultrospec 2100 Pro



Amersham Ultrospec 2100 Pro UV/Vis spectrophotometer, 190-900 nm, 8-position cell changer, on board LCD microprocessor display

Rating: Not Rated Yet

**Price**

Sales price \$4,250.00

[Ask a question about this product](#)

Manufacturer [Amersham](#)

### Description

This Amersham Ultrospec 2100 Pro UV/Vis spectrophotometer has a 190-900 nm range, 8-position cell changer, on board LCD microprocessor display, automatic wavelength calibration on start up, serial and parallel Centronics interfaces, & provides stored routines for DNA, RNA and oligonucleotide quantitation and purity check.

Wavelength scan with zoom, absorbance changes with time, reaction rate determinations, and standard curves can be displayed and printed out. User defined equations can be entered using the multi-wavelength mode and up to 18 methods can be saved. Longer lasting Xenon lamp configuration.

Nice compact design & beautiful yellow color.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)