

## Arrow JR 4000 Overhead Stirrer



Light duty; 4000 rpm max

Rating: Not Rated Yet

**Price**

Sales price \$300.00

[Ask a question about this product](#)

### Description

Reliable for light-duty mixing applications. Robust construction, use-friendly design. Batch sizes up to 20 liters. With sturdy mounting rod and grounded plug and power cord for safety.

### Key Features and Specifications:

Motor Power: 1/25 horsepower, delivering consistent performance with a maximum speed of 4,000 RPM.

Torque: Rated at .66 in-lbs, providing the necessary force for effective mixing.

Propeller Shaft: 3/8" diameter by 12" long, made from durable stainless steel (S-316-12).

Propeller: 2 1/2" diameter (#25B), designed for efficient mixing.

Mounting Rod: 9/16" diameter by 10" long (#15) for stable and secure setup.

Power Supply: 120 VAC 60 Hz, with a 6ft power cord and grounded plug for safety and convenience.

Gear Ratio: Direct, ensuring straightforward and reliable operation.

## Shakers/Stirrers/Hotplates: Arrow JR 4000 Overhead Stirrer

---

Shipping Weight: 6 lbs (3.2 kg)

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)