

## Barnstead Thermolyne Type 45600 Cellgro Stirrer





With Soft Start/Soft Stop speed control with a gradual change of 5.5 rpm.

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Manufacturer [Barnstead](#)

#### Description

The Barnstead Thermolyne Type 45600 Cellgro Stirrer, Model # S45625, is a five-position magnetic stirrer designed for precise, controlled stirring to optimize cell growth in laboratory settings. Featuring a Soft Start/Soft Stop speed control with a gradual change of 5.5 rpm to minimize cell trauma, it can accommodate up to five 1000 ml beakers. Its robust design supports operation in a wide ambient temperature range of -20°C to 60°C, making it suitable for benchtop or incubator use, with an auxiliary receptacle for operating multiple units from one outlet.

#### Specifications:

Model: S45625

Stirring Positions: 5

Capacity: Up to 1000 ml per position

Speed Control: Soft Start/Soft Stop (5.5 rpm gradual change)

Ambient Temperature Range: -20°C to 60°C

Dimensions: 7 x 7 x 9.5 in. (inside), 8 x 12 x 13 in. (outside)

Features: Auxiliary receptacle for multi-unit operation

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)