

Beckman 3.7 Swing Bucket Rotor



The Beckman GH 3.7 Swing Bucket Rotor is an unshielded, four-place rotor with stainless-steel rotor yoke and removable aluminum swinging buckets.

Rating: Not Rated Yet

Price

Sales price \$500.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

Description

The **Beckman GH 3.7 Swing Bucket Rotor** is an unshielded, four-place rotor with stainless-steel rotor yoke and removable aluminum swinging buckets. The buckets are interchangeable with MicroPlus multiwell plate carriers (PN # 362394) for spinning microtiter plates or microfuge tubes using rack inserts. Mostly used to rapidly sediment protein precipitates, large particles, cells and cell debris, for binding studies and for separating serum from whole blood.

[Biotech Equipment Sales, Inc.](#)
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com