

Beckman Allegra X30R Centrifuge



Centrifuges: Beckman Allegra X30R Centrifuge

Benchtop; w 4-place swinging bucket rotor SX 4400

Rating: Not Rated Yet

Price

Sales price \$5,250.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

Description

The X-30R refrigerated model cools from room temperature to 4° C in fewer than ten minutes and has an operating temperature range of 2° C to 40° C (35° F to 104° F). The X-30R provides excellent performance and easy operation in a compact design that saves valuable time and laboratory space with the highest speed and g-force in its class (1.6L max capacity, 18,000 RPM, 29,756 x g). Simple and quick run setup. Versatile, with an extensive rotor library for multiple applications, plus a compact Design of only 18" wide (46cm). Rotor: 4-place swinging bucket rotor SX 4400

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com