

Beckman Coulter Allegra 25R



NEW Beckman Coulter Allegra 25R high-speed refrigerated benchtop centrifuge

Rating: Not Rated Yet

Price

Sales price \$6,350.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

!!!! - RECENTLY SOLD - !!!!

This Beckman Coulter Allegra 25R - BRAND NEW IN CRATE - refrigerated large capacity refrigerated benchtop centrifuge is microprocessor controlled and can get up to 15,000 rpm, has a temp range of -2° to 40°C , can hold up to 4 x 500 ml (at 5700 rpm) w/the large capacity swinging bucket rotor, 208V. Rotor not included & can be ordered from Beckman Coulter.

Please call us at (650) 871-5707 to check availability of this product!