

## Beckman Coulter Allegra V-15R Centrifuge



Allegra V-15R refrigerated centrifuge microplate package, includes instrument and VS 4.750 96 swinging-bucket rotor.

Rating: Not Rated Yet

**Price**

Sales price \$10,950.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

Model: Allegra V-15R

120VAC, 16A, 60Hz

Rotor VS 4.750

Features

- 
- ? 56 dBA at full speed
- Intuitive, uncomplicated panel
- Bright, at-a-glance display
- Automatic rotor recognition

- Rotor life tracking
- Built-in imbalance sensors
- BioSafe rotors
- 50 programmable runs
- 10 rotor configurations
- Comprehensive selection of adapters for swinging-bucket and fixed-angle rotors
- Speeds up to 13,500
- Spins up to 4 x 750 mL
- Force up to 20,412 x g
- Up to 28% smaller than comparable models
- 60.5 centimeters wide and 63.5 centimeters deep

### Specifications

Timed Run To 99 hrs: 59 mins or Continuous run (hold)

Max Capacity Volume Range: 0 – 3 L

Temperature Control: -10° C – +40° C

Temperature Range(s): From -10°C to +40°C

Display: Backlit LCD

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)