

Beckman Coulter Avanti J25i



Beckman Coulter Avanti J25i Superspeed floor model centrifuge

Rating: Not Rated Yet

Price

Sales price \$9,750.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

This Beckman Coulter Avanti J25i Superspeed high performance floor model centrifuge, microprocessor controlled to 25,000 rpm, programmable, -20C to 40C, induction drive, not including rotors

!!!! - RECENTLY SOLD - !!!!