


## Beckman Coulter Avanti JLA-8.1000 Rotor, Manufactured 2022





Max 8000 RPM

Rating: Not Rated Yet

**Price**

Sales price \$9,750.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

General purpose, light weight rotor for large-volume processing, as well as pelleting of bacteria, subcellular organelles, viruses and precipitates.

Specifications

JLA-8.1000

Part no.: 3636885

No. tubes x volume (ml): 6x1,000

Max speed: 8,000RPM

Max force: 15,970x g

K factor: 2,482

Instrument used with: JXN-26, J-26S, J-HC

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)