

Beckman Coulter DU530 Spectrophotometer



The DU-530 is a simple but great spectrometer that comes with various built-in programs that are easy to use.

Rating: Not Rated Yet

Price

Sales price \$3,250.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

The DU-530 is a simple but great spectrometer that comes with various built-in programs that are easy to use. This system uses a single cell module that uses a wide variety of cell holders such as, the 1cm standard, 16mm test tube, 5cm and 10cm rectangular and a 50µl microcell.

Built-in Software Options Include -

- * Fixed Wavelength Mode
- * Scanning Mode
- * Single Component Analysis
- * Time Based Kinetics
- * Protein Analysis
- * Nucleic Acid Analysis with DNA/RNA/Oligo Quantitation & Purity Checks

Specifications -

Cell Module	1cm
Wavelength Range	190 - 1100 nm
Programs for Methods	100 user data storage
Interface	RS-232/ Parallel Port
Software	Multi Language

