

## Beckman Coulter JA-17 Fixed-Angle Rotor



Rotor Accommodates: 14 50 mL Tubes; Maximum Speed: 17,000 g

Rating: Not Rated Yet

**Price**

Sales price \$2,750.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

### Description

Aluminum, 14-place, and fixed-angle rotor is for use in the Avanti J-25 XP and XPIBK393124 series and J-EBK369001 and J6-MIBK360291 centrifuges. Ideal for high-force, large-volume applications, including pelleting bacteria, harvesting cell membranes, and separating subcellular organelles.

Applications: Harvesting bacteria and cell membranes, processing tissue homogenates, and separating cell particulates

Rotor Accommodates: 14 50 mL Tubes

Maximum Speed: 17,000 g

Roter Capacity: 700 mL

Approximant Acceleration Time: 2:30

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)