

## Beckman Coulter JA-20.1 Fixed-Angle Rotor



32-place fixed-angle rotor

Rating: Not Rated Yet

**Price**

Sales price \$1,450.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

### Description

This aluminum 32-place, fixed-angle rotor is for use with Avanti J-26 XP and XPIBK393124 series and J-EBK369001 and J6-MIBK360291 centrifuges.

Ideal for high-force, large-volume separation of particles with 100S or larger sedimentation coefficients.

Application: High-Force, Large-Volume Separation of Particles with 100S or larger sedimentation coefficients.

Material: Aluminum

32-Place Roter

Fixed-Angle Rotor Accommodates up to 32 15 mL Tubes

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)