

Beckman Coulter JLA-10.500, Never Used



6 x 600 mL; 10,000 rpm; 18,600 x g

Rating: Not Rated Yet

Price

Sales price \$7,500.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

Model: JLA-10.500

Fixed-Angle Rotor, Aluminum

Includes 6 Biosafety Liquid Canisters

Major applications: Large volume for initial processing of bacteria or other cells from fermenters, clearing cell debris from homogenates, or pelleting subcellular organelles and protein precipitates. For use in Avanti J Series and J2 Series centrifuges.

6 x 600 mL

10,000 rpm

18,600 x g

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office

650-276-7487 fax
sales@biotechequipmentsales.com