

## Beckman Coulter JS 4,2 Swinging Bucket Rotor



Windshielded, six-place rotors with aluminum rotor yoke and removable aluminum swinging buckets for use in the Avanti™ J6-MI centrifuge. Ideal for rapid sedimentation of protein precipitates, large particles, cells, and cell debris. Suitable fo

Rating: Not Rated Yet

**Price**

Sales price \$4,500.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

### Description

Windshielded, six-place rotors with aluminum rotor yoke and removable aluminum swinging buckets for use in the Avanti™ J-HC and J6 Series centrifuges. Ideal for rapid sedimentation of protein precipitates, large particles, cells, and cell debris. Suitable for binding studies and separating serum from whole blood.

This rotor has a capacity of 6L or six blood bags. The buckets are interchangeable with microplate carriers (sold separately); each carrier holds three microtiter plates.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)