

Beckman Coulter MLA-80 Fixed-Angle Rotor



MLA-80 Rotor Package, Fixed Angle Rotor (26 Angle), Titanium.

Rating: Not Rated Yet

Price

Sales price \$5,450.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

Major applications: Large volume differential sedimentation (pelleting) of subcellular organelles and viruses. For use in Optima MAX-XP Benchtop Ultracentrifuge and Optima MAX High-Capacity Personal Micro-Ultracentrifuge.

Specifications

- Platform: Tabletop UCF
- Rotor Type: Fixed-Angle (Ultracentrifugation)
- Maximum RPM: 80,000 rpm

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com