

## Beckman Coulter Optima L-90K



Beckman Coulter Optima L-90K refrigerated floor model ultracentrifuge to 90,000 rpm. Very low hours!!

Rating: Not Rated Yet

**Price**

Sales price \$22,500.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

Description

**!!!SOLD!!!**

This Beckman Coulter Optima L-90K refrigerated floor model ultracentrifuge is capable of generating 694,000 x g at speeds up to 90,000 rpm, the Optima L-90 K Ultracentrifuge enables you to perform more separations in less time. This versatile floor model operates with a broad range of superb rotors, including zonal and continuous flow for large-volume separations. The L-90 K offers the reliability of a vacuum-encased induction drive, the simplicity of user-friendly, microprocessor-based control, and environmentally-friendly cooling systems that eliminate the need for CFCs and other harmful liquid refrigerants.

Ours only has 600 million revolutions on it (about 300 hours total use) and has been well-maintained. Has teh original logbook. We also have some barely used ultra rotors listed in our rotor section. Comes with our 90-day warranty.

**Recently Sold: Beckman Coulter Optima L-90K**

---

Service contract available in the Bay Area and many other regions.

**Please call us at (650) 871-5707 to check availability of this product!**