

## Beckman Coulter Type 70 Ti Fixed-Angle Titanium Rotor - NEW IN BOX!



Maximum RPM: 70,000 rpm; Max g-Force: 504,000.00

Rating: Not Rated Yet

**Price**

Sales price \$11,950.00

[Ask a question about this product](#)

Manufacturer [Beckman Coulter](#)

#### Description

Product No: 337922

For use in instruments classified: HRS; Major Applications: Differential centrifugation of subcellular particles.

#### Specifications

Platform: Floor Prep

Rotor Type: Fixed-Angle (Ultracentrifugation)

Angle: 23

Maximum RPM: 70,000 rpm

Max g-Force: 504,000.00

k-Factor: 44

Number of Tubes: 8x39 mL / 25x89 mm(8x39 mL / 1x4 in)

Caps and Spacers: 342699 ? 343448 ? 344635 ? 361669

Nominal Capacity: 312 mL

rMAX: 91.9

rMIN: 39.5

rAV: 65.7 mm

BioSafe: Bio-Enhanced

Materials: Titanium

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)