

## Beckman GS-6R



Beckman GS-6R refrigerated benchtop centrifuge w/GH 3.8 4-place swinging bucket rotor

Rating: Not Rated Yet

**Price**

Sales price \$2,875.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

### Description

This Beckman GS-6R refrigerated benchtop centrifuge w/GH 3.8 4-place swinging bucket rotor has a capacity of 4 x 750 ml, max speed w/rotor of 3750 rpm, temp range of -5 to 25C, & 30 minute timer. Many inserts in stock. Multiple units are available.