

Beckman JA-20 Fixed Angle Rotor



Major applications: Harvesting bacteria and cell membranes, processing tissue homogenates, separating cell particulates.

Rating: Not Rated Yet

Price

Sales price \$1,500.00

[Ask a question about this product](#)

Manufacturer [Beckman](#)

Description

- Item JA-20 Rotor, Fixed Angle, Aluminum, Biosafety Lid, 8 x 50 mL, 20,000 rpm, 48,400 x g
- Company [Beckman Coulter](#)
- Price Inquire
- Catalog Number 334831
- Quantity Each

Specifications/Features

Association:

JA-20

Fixed-Angle:

Fixed Angle

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com