



## Bio-Rad GelDoc XR+ Gel Photo Documentation System



For fast and easy quantitation of a variety of nucleic acid gels, protein gels, and blots.

Rating: Not Rated Yet

**Price**

Sales price \$8,450.00

[Ask a question about this product](#)

Manufacturer [BioRad](#)

Description

The GelDoc XR+ is a high resolution gel documentation system that allows for fast and easy quantitation of a variety of nucleic acid gels, protein

gels, and blots.

Includes:

USB CCD camera and a motorized zoom lens with software driven autofocus,

Universal Hood II darkroom with built-in white light epi-illumination and UV transillumination

2-position manual filter changer with standard Amber/ethidium bromide filter (548-630nm)

Optional white light transillumination conversion screen

Windows 10 computer with software installed, monitor, keyboard and mouse

Power req: 100-230V, 50/60Hz

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)