

## BioTek EL406 Combination Washer Dispenser 406PSUB3



Fast, full (96-, 384-well) plate-washing along with three reagent dispensers in one compact instrument.

Rating: Not Rated Yet

**Price**

Sales price \$25,750.00

[Ask a question about this product](#)

Manufacturer [BioTek](#)

### Description

Capabilities include washing and dispensing for standard assays, plus biomagnetic separation and vacuum filtration protocols. The Dual-Action™ manifold optimizes washing for loosely adherent cell monolayers, while the patented Ultrasonic Advantage™ helps prevent blockages in the manifold tubes before washing, to reduce assay failure. Syringe pump and peristaltic pump dispensers each offer unique advantages, and the EL406 offers both, along with superior plate washing, on a single platform.

- 96-, 384- well plate compatible
- Gentle processing for cell-based assays
- Bead washing: magnetic and polystyrene

Ask about optional plate-stacker.

[Biotech Equipment Sales, Inc.](#)  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)