

## BioTek ELx405 UVS Select Microplate Washer



96- and 384-well; with Ultrasonic Advantage.

Rating: Not Rated Yet

### Price

Sales price \$4,500.00

[Ask a question about this product](#)

Manufacturer [BioTek](#)

### Description

With wash buffer switch and Ultrasonic Advantage. Ask about bottles, vacuum pumps, etc. to complete your system.

#### Specifications:

- Microplate Types: 96 and 384-well, Low profile and standard height
- Onboard Software: Create, edit and run multiple protocols
- Soaking: Programmable up to 60 minutes
- Manifold Types: Two manifolds, one with 96 aspirate tubes and the other with 96 dispense tubes
- Volume Range: 50 to 3000 microliter increments
- Wash Cycles: 1 -10
- Buffer/reagent Selection: Auto switching module for up to 4 buffers (option)
- Dispense Precision: 24% CV: 300 µL/well (96-well washing)
- Residual Volume:
- Wash Speed: 96-wells, 300 µL/well, 3 cycles: 384-wells, 100 µL/well, 3 cycles
- Flow Rates: High flow to low flow
- Dimensions (D x W x H): 17" x 14" x 10" (43.2 x 35.6 x 25.4 cm)
- Weight: 30 lb. (13.5 kg)

---

Biotech Equipment Sales, Inc. based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

---

[Biotech Equipment Sales, Inc.](#)  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)