

Blue M Oven



4.8 cu ft.; temp limit 260C

Rating: Not Rated Yet

Price

Sales price \$1,625.00

[Ask a question about this product](#)

Manufacturer [Blue M](#)

Description

Designed for laboratory applications such as moisture drying, curing, baking, aging and general testing. The ovens have been designed for convenience and appearance in the workplace and include ergonomic handles and controls, with a powder-coated steel exterior. Components are high quality, to ensure years of reliable service.

Mechanical convection; horizontal airflow; stainless steel interior. High-quality digital PID temp controller; setpoint operation; 2 shelves.

Specifications

Ovens/Furnaces: Blue M Oven

Model No.: LO - 136 - P

Interior dimensions (WxDxH): 550 x 450 x 550 mm / 21.6 x 17.7 x 21.6 in.

Exterior dimensions (WxDxH): 750 x 690 x 1045 mm / 29.5 x 27.2 x 41.1 in

Workspace: 4.8 cu ft

Temp range: Ambient + 10C: + 260C

Temp uniformity: ? + 3C

Control tolerance: ? +1C

Heating rate: ?100 min (ambient temp to max temp)

Power input: 230V / 1PH / 60Hz

FLA: 9.5A

Rated heat input: 2.0KW

SCCR: 5KVA

NEMA 1 Enclosure

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com