

Boekel Scientific MicroCooler II Model 260010



Temp Range: 18°C below ambient

Rating: Not Rated Yet

Price

Sales price \$550.00

[Ask a question about this product](#)

Manufacturer [Boekel](#)

Description

Ideal for maintaining the 14°C used in ligation reactions or 17°C for storage applications, as well as maintaining samples at ice bucket temperature and cooling blood samples prior to coagulation testing. Compact design conserves bench space. The watertight well can function as a small refrigerated bath and features a "press to set" microprocessor controller.

Specifications:

Temp Range: 18°C below ambient

Accuracy: +/- 1°C of set point

Capacity: 1 rack

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com

