

## Buchi R-210 with Specialized B-495 Bath



R-210 Rotary Evaporator with Specialized B-495 self replenishing bath.

Rating: Not Rated Yet

**Price**

Sales price \$3,450.00

[Ask a question about this product](#)

Manufacturer [Buchi](#)

### Description

**!!!! - RECENTLY SOLD - !!!!**

- Compact assembly with small footprint on your lab-table
- Easy handling with the electrical quick-action jack activated by simply pressing a button
- Sealing system: A continuous vapor duct prevents cross-contamination and is particularly easy to dismantle and clean. The vacuum seal comprises a PTFE composite material and provides long service life and effective sealing. The proven combi-clip elegantly solves the problem of simple flask removal and the dismantling of the vapor duct
- The heating bath covers a wide temperature range from 20-100°C. A large, clear graphic display continuously shows the actual temperature and the temperature setpoint. Its small volume ensures that warm-up times are very quick. Thanks to the wireless power supply, it is not necessary to make or break connections when handling the bath. The high alloy steel pan is easy to clean. An optional protective shield can easily be fitted to the edge of the bath.
- If required, all glass parts (except the evaporating flask) can be supplied with a special plastic coating. The P+G coating PLASTIC+GLASS\* protects against mechanical damage and provides efficient protection against glass splinters in the unlikely event of a breakage. (in option )
- Total flexibility with his modular approach

\* PLASTIC+GLASS ( P+G) is a unique protective layer for glass components. It offers improved mechanical damage resistance and increases protection against broken glass. It also makes sure that the solvent in the receiving flask is not spilling, if the glass is damaged.

With a heating capacity of 1300 W, the B-495 provides enough energy for a fast distillation. In addition the B-495 can be used as a water bath up to 100°C. The B-495 is additionally equipped with a automatic water level control (bath replenishment).

Recently Sold: Buchi R-210 with Specialized B-495 Bath

---

Complete with dry ice condenser and two 1 Liter flasks.

-!!! SOLD !!!-