

Buchi R-220 Pro Rotary Evaporator



With 20 L evaporating flask, two receiving flasks, single condenser.

Rating: Not Rated Yet

Price

Sales price \$33,500.00

[Ask a question about this product](#)

Manufacturer [Buchi](#)

Description

Scale-up your laboratory process with the Rotavary Evaporator R-220 Pro.

The key with this system is that it has the cart below it that has both the chiller and the vacuum pump, making it a "plug and play" system. Includes the Buchi F-325 Recirculating Chiller and the Buchi V-600 Vacuum Pump. With 20 L evaporating flask, two receiving flasks, and a single condenser.

Evaporators/Concentrators: Buchi R-220 Pro Rotary Evaporator

Distillation rate: up to 12 L Ethanol/hour

Heating bath temperature: up to 180 °C

Sample size: max. 10 L per batch

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com