

Caliper Life Sciences TurboVap LV Concentration Workstation



Automated alternative to rotary and vacuum-style evaporators.

Rating: Not Rated Yet

Price

Sales price \$3,750.00

[Ask a question about this product](#)

Description

Designed to work with tube sizes from 1.5mL to 35mL and evaporate/concentrate up to 50 samples simultaneously. Rack holds 15 mL tubes.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com