

Cole-Parmer Digital Gear Pump 75211-30



The IP23-rated enclosure protects internal components from damage. Digital optical encoder feedback signal keeps speed control to $\pm 0.3\%$.

Rating: Not Rated Yet

Price

Sales price \$700.00

[Ask a question about this product](#)

Manufacturer [Cole-Parmer](#)

Description

The IP23-rated enclosure protects internal components from damage. Digital optical encoder feedback signal keeps speed control to $\pm 0.3\%$.

Pump head not included, can be ordered from manufacturer.

Specifications

Temp range	104°F (40°C) maximum
Duty cycle	Continuous
Motor type	Permanent magnet DC (PMDC), IP23
Power (Amps)	3
Power (Hz)	50/60
Dimensions	7-3/4"W x 7-1/4"H x 11-1/2"D
Motor (hp)	1/10
Motor rpm	60 to 3600
Power (VAC)	115
Brand	Cole-Parmer
Manufacturer number	75211-30

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com