

Edvotek EdvoCycler



The EdvoCycler™ is a standalone classroom PCR machine that is easy to use! The 0.2 ml tube block has room for up to 25 student samples and comes preprogrammed with all Edvotek PCR kit protocols. These programs may be modified or deleted, plus there's extra

Rating: Not Rated Yet

Price

Sales price \$850.00

[Ask a question about this product](#)

Manufacturer [Edvotek](#)

Description

The EdvoCycler™ is a standalone classroom PCR machine that is easy to use! The 0.2 ml tube block has room for up to 25 student samples and comes preprogrammed with all Edvotek PCR kit protocols. These programs may be modified or deleted, plus there's extra memory for more. The vivid 7-line LCD displays all program parameters simultaneously on a single screen. A heated oil-free lid makes operation a snap.

Holds 25 x 0.2 ml PCR samples

- Heated oil-free lid with magnetic latch
- Pre-programmed Edvotek PCR protocols
- Vivid 7-line LCD display with live program information
- Standalone machine - no PC required!
- Temperature Range: 4°C - 99°C
- Maximum Ramp Rate: 3°C/sec.
- Dimensions: 16 x 8.5 x 7"

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com