

## Eppendorf 5415D Centrifuge



The 5415 D has all the tried and tested technical features, including the almost indestructible maintenance-free drive. It goes without saying that the easy-to-follow digital display shows all significant parameter settings at a glance.

Rating: Not Rated Yet

**Price**

Sales price \$800.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

### Description

The 5415 D has all the tried and tested technical features, including the almost indestructible maintenance-free drive. It goes without saying that the easy-to-follow digital display shows all significant parameter settings at a glance.

The centrifuge can be operated using the same 24-position standard rotor for 1.5 ml and 2 ml tubes or with the 36-position special rotor for 0.5 ml tubes. With or without refrigeration - the right choice is always 5415!

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**