

### Eppendorf 5702RH Benchtop Refrigerated/Heated Centrifuge



Includes A-4-38 swinging-bucket rotor with 15 mL and 50 mL conical-tube adapters.

Rating: Not Rated Yet

**Price**

Sales price \$5,500.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

**Description**

Eppendorf 5702RH is a compact, refrigerated/heated benchtop centrifuge designed for routine cell culture, blood, and general laboratory centrifugation. With its A-4-38 swing-bucket rotor, it accommodates tubes from about 1.5–100 mL; the included adapters allow use of 15 mL and 50 mL conical tubes.

**Key Features/Specifications**

Maximum speed: 4,400 rpm

Maximum RCF: 3,000 × g

Swinging-bucket capacity: 4 × 100 mL with A-4-38 rotor

Temperature control: refrigerated/heated, approximately –9°C to +42°C, allowing temperature-sensitive samples to be processed at controlled temperatures.

Tube versatility: supplied configuration includes adapters for 15 mL and 50 mL conical tubes, with the rotor supporting a broad range of tubes from approximately 1.5–100 mL.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**

**226 Miller Ave.**

**South San Francisco, CA 94080**

**650-871-5707 office**

**650-276-7487 fax**

[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)