

Eppendorf 5425 Centrifuge



Centrifuges: Eppendorf 5425 Centrifuge

Includes Rotor FA-24x2

Rating: Not Rated Yet

Price

Sales price \$1,950.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

Whisper quiet 24-place standard microcentrifuge with outstanding feature set for broad range of applications.

With touchpad control panel for easy cleaning.

Includes Rotor FA-24x2

Non-refrigerated

120 V, 50 – 60 Hz

Max speed (RCF): 21,300 × g

Max RPM: 15,060

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com