

## Eppendorf F45-24-11 Fixed Angle Rotor



Eppendorf F45-24-11 Fixed-Angle Rotor

Rating: Not Rated Yet

**Price**

Sales price \$325.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

### Description

Eppendorf F45-24-11 Fixed-Angle Rotor for use with Eppendorf Centrifuges 5702, 5702R, 5702RH, 5415D, 5415R.

### Specifications -

- .Capacity: 24 x 1.5/2.0 mL tubes
- .Max. Rpm: 13,200 rpm
- .Max. rcf: 16,110 x g
- .Autoclavable (121°C, 20 minutes)

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)