

Eppendorf F45-30-11 Fixed Angle Rotor



High capacity rotor for 30 x 1.5/2.0 mL tubes Max. rcf: 20.800 x g (14.000 rpm) Rotors, rotor lids and adapters are autoclavable (20 min, 121°C) Adapters available for - 0.2 mL PCR tubes - 0.4 mL centrifuge tubes - 0.5 mL mic

Rating: Not Rated Yet

Price

Sales price \$350.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

- High capacity rotor for 30 x 1.5/2.0 mL tubes
- Max. rcf: 20.800 x g (14.000 rpm)
- Rotors, rotor lids and adapters are autoclavable (20 min, 121°C)
- Adapters available for
 - 0.2 mL PCR tubes
 - 0.4 mL centrifuge tubes
 - 0.5 mL micro test tubes
 - 0.6 mL Microtainer® tubes

[Biotech Equipment Sales, Inc.](#)
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com