

Eppendorf Mastercycler Personal



Smaller 25-place thermal cycler, 4 to 99C w/simple programming and LCD display.

Rating: Not Rated Yet

Price

Sales price \$1,850.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

!!!SOLD!!!

Mastercycler Personal is a high quality, compact thermal cycler. Along with speed, a heated lid, and precise block uniformity, its universal block with 25-well format can accommodate 25 x 0.2 mL tubes, 16 x 0.5 mL tubes or microplates in a 5 x 5 grid format. Programming includes adjustable ramps, time and temperature increments, incubation mode, link function, and pause function.