

Eppendorf / New Brunswick Innova 42R Refrigerated Incubator Shaker



Shaking speed range from 25-400 rpm, temperature control from 20°C below ambient to 80°C, 0.1-99.9 hour timer

Rating: Not Rated Yet

Price

Sales price \$7,500.00

[Ask a question about this product](#)

Manufacturer [Eppendorf](#)

Description

Mfg no. M1335-0014, Shaking speed range from 25-400 rpm, temperature control from 20°C below ambient to 80°C, 0.1-99.9 hour timer, features triple-eccentric counterbalanced drive, solid state DC brushless motor, 4 storable programs with 15 programmable steps each, 1in. orbital throw, 18x18in. shaker platform dimensions, comes with optional culture drawer, 120V, 60Hz, 12.5A.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com