

## Fisher Scientific Electrophoresis Power Supply FB300



Maximum voltage: 300V

Rating: Not Rated Yet

**Price**

Sales price \$325.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

Features

Running up to three gels at once at either constant current or constant voltage  
Stackable, solid, compact power source  
Integral timer and auto shutoff prevents fear of losing valuable information  
Ideal for mini SDS-PAGE, DGGE/CDGE and DNA/RNA submarine electrophoresis  
Maximum voltage: 300V  
Maximum current: 400mA  
Maximum power: 75w

Three jacks

Internal timer: 0 to 999 minutes

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)