

Fisher Scientific Isotemp 100L Gravity Incubator



100-liter capacity; temperature range from ambient +5°C to 75°C

Rating: Not Rated Yet

Price

Sales price \$1,750.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

The Fisher Scientific Isotemp 100L Gravity Incubator (Cat. No. 15130514) is a reliable and efficient solution for microbiological and laboratory applications, offering precise temperature control and uniformity through gravity convection technology. With a 100-liter capacity and a temperature range from ambient +5°C to 75°C, it supports a variety of routine lab tasks. Its stainless-steel construction, digital display, and internal glass door ensure durability, ease of use, and sample visibility without heat loss, making it ideal for research and testing environments.

Specifications:

Model: Isotemp 100L Gravity Convection Incubator (Cat. No. 15130514)

Capacity: 100 L (3.53 cu. ft.)

Temperature Range: Ambient +5°C to 75°C

Convection Technology: Gravity convection

Exterior Material: Stainless steel (AISI 430 / 1.4016)

Interior Features: Internal glass door for sample viewing

Control Type: Digital with fluorescent display

Shelving: Flexible shelving system, 2 included; 55 lbs shelf weight capacity

Additional Features: Access port with silicone plug, stackable design, intuitive controls

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com