

Fisherbrand Isotemp FBG4505GA Laboratory Refrigerator - from 2024!



Two sliding doors; 8 adjustable shelves.

Rating: Not Rated Yet

Price

Sales price \$4,150.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

The Fisherbrand Isotemp FBG4505GA is an ENERGY STAR certified upright laboratory refrigerator designed for clinical and research applications in the United States and Canada. This model features two sliding glass doors and eight adjustable shelves, offering a spacious 45.8 cu. ft. (1297L) capacity. With a quiet operation below 53 dBA and sustainable R290 refrigerant, it ensures reliable temperature control and eco-friendly performance.

Specifications:

Capacity: 45.8 cu. ft. (1297L)

Refrigerators: Fisherbrand Isotemp FBG4505GA Laboratory Refrigerator - from 2024!

Shelves: 8 adjustable
Temperature Range: 3°C to 7°C
Refrigerant: R290 (GWP: 3)
Defrost: Automatic (heat-free)
Energy Consumption: 2.41 kWh/day
Dimensions (D x W x H) Exterior: 36.6 x 56.5 x 78.5 in. (93.1 x 143.5 x 199.4 cm)
Dimensions (D x W x H) Interior: 26.0 x 52.5 x 58.0 in. (66.0 x 133.4 x 147.3 cm)
Electrical Requirements: 115V, 60Hz, NEMA 5-15 plug
Certifications: cULus, ENERGY STAR

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com