

Fisherbrand Isotemp GDP 10 Water Bath



10L Deluxe General Purpose Water Bath

Rating: Not Rated Yet

Price

Sales price \$450.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

Model FSGPD10

Ideal for industrial, clinical pharmaceutical or biomedical applications.

Features:

- Smaller footprint, compared to previous models
- Advanced temperature controller is designed for extended functionality
- New audible alarms
- Conveniently save commonly used settings with temperature presets
- Simple to operate and monitor with icon-based graphical display
- Enhanced hinge design supports open lid configuration and easy lid removal
- Help prevent bath damage and overheating with low-fluid protection
- Compatible with Thermal beads (sold separately)

Specifications:

- Chamber Material: Stainless Steel
- Description: 10L Deluxe General Purpose Water Bath
- Display: Monochrome LCD
- Includes: Polycarbonate Gable Cover, Diffuser Tray
- Model: Isotemp GDP 10
- Temperature Stability: $\pm 0.2^{\circ}\text{C}$
- Dimensions (D x W x H): Interior 11.9 x 13 x 5.99 in. (301 x 330 x 150 mm)
- Amperage: 6.0-7.0 / 2.8-3.4 (120V/230V)
- Controller Type: Digital

Baths/Circulators: Fisherbrand Isotemp GDP 10 Water Bath

Capacity (Metric): 10 L
Temperature Range: Ambient to 100°C
Heating Capacity: 800 w
Material: Epoxy Powder-coated Steel
Temperature Presets: 3
Temperature Uniformity: $\pm 0.2^{\circ}\text{C}$
Dimensions (L x W x H): 15.5 x 15.1 x 9.2 in. (393 x 383 x 233 mm)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com