

Fisherbrand Model 505 Sonic Dismembrator



Hertz: 20 kHz (Ultrasonic Frequency); Volume Range (Metric): 0.5 to 500 mL

Rating: Not Rated Yet

Price

Sales price \$3,750.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

Description

High power, programmable system ideal for nanotechnology, cell disruption and homogenization applications.

Specifications

Hertz: 20 kHz (Ultrasonic Frequency)
Volume Range (Metric): 0.5 to 500 mL
Height (English): 9 in.

Shakers/Stirrers/Hotplates: Fisherbrand Model 505 Sonic Dismembrator

Width (English): 18 in.

Timer: 1 sec. to 10 hrs.

Wattage: 500 W

Dimensions (L x W x H): 15 x 18 x 9 in.

Length (English): 15 in.

For Use With (Application): Ideal for nano particle dispersion, mixing and homogenization

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com