

## Fisherbrand Nutating Mixer - Variable Speed



## Shakers/Stirrers/Hotplates: Fisherbrand Nutating Mixer - Variable Speed

---



Usable in non-CO2 incubators, refrigerators and chambers.

Rating: Not Rated Yet

**Price**

Sales price \$450.00

[Ask a question about this product](#)

Manufacturer [Fisher Scientific](#)

**Description**

Delivers gentle uniform mixing of fluids in tubes, plates, culture flasks and dishes. (Includes rubber mat, not pictured.)

Model 88861043

Variable speed from 2 to 60rpm offers gentle continuous mixing of samples. Low-profile, compact design enables placement at every workstation. 12 x 10.5 in. platform accommodates a variety of tubes and vessels. Display will show last used settings, even after power has been turned off.. Operating temperature range - 5° to 40°C, usable in non-CO2 incubators, refrigerators and chambers. Continue to run by last set speed and the remaining time when the power recovered. Continuous or timed operation with automatic switch-off. Audible and visual alarm: In timed mode, when the time reaches zero, alarm will sound five times, then the time window shows "End" and flickers all the time. Alarm can be stopped by touch any key

**Specifications**

Rocking Angle: 20°

Speed Control Accuracy: ±1 rpm

Electrical Requirements: 100/240 V, 1 A, 9.5 W, 50/60 Hz

Temperature Range (Metric): 5°C to 40°C, usable in non-CO2 incubators, refrigerators and chambers

Load Bearing Capacity (English): 1.8 lb.

Plug Type: US Plug

Material: Aluminum Platform

Dimensions (L x W x H) 16.3 x 11.5 x 10.7 in. (415 x 293 x 273 mm)

Description: 20° Tilt Nutating Mixer - Variable Speed

Includes: Dimpled rubber mat and six tube securing bands

Load Bearing Capacity (Metric): 0.8 kg

Speed Range: 2 to 60 rpm

Shipping Weight (English): 13.2 lb.

Alarms: Audible and Visual

Display: LED

Timer Range: 0 min. to 99:59 hr.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)