

## Forma 416 Shaker



Forma 416 large-capacity benchtop orbital shaker

Rating: Not Rated Yet

**Price**

Sales price \$1,750.00

[Ask a question about this product](#)

Manufacturer [Forma](#)

Description

**!!!SOLD!!!**

This Forma 416 large-capacity benchtop orbital shaker has an 18" x 24" universal platform which can accept many different flask clamps. Microprocessor control of speed from 25 to 525 rpm, and of time to 200 hrs.

**Please call us at (650) 871-5707 to check availability of this product!**